

Large Capacity Suction System



Features

- ▶ **All in One, saving 50% space**
Lafil 200 suction system integrates vacuum source, waste bottle and suction kit into one single system. Saving 50% bench space.
- ▶ **Large volume, effort-saving**
Equipped with 4L large-capacity PP bottle, saves efforts of repetitive waste draining. (Optional 3.5L glass bottles are also available)
- ▶ **Safety features**
Multiple safety featured designs protect both the operator and the machine.
(1) Easier and safer waste drain with quick-coupling bottle connectors.
(2) Tripping-proof and spill-proof fence-like bottle platform.
(3) Overflow protection prevents liquid from entering the pump
- ▶ **Eco control feature (Lafil 200^{eco})**
Lafil 200^{eco}'s pump pauses when maximum vacuum is detected and resumes operation to maintain a constant and stable suction. Making Lafil 200^{eco} more energy-saving and quieter machine to work with.
- ▶ **International certification**
CE certification
IP 41 protection
RoHS certification
- ▶ **One-Year Warranty (Lafil 200 / Lafil 200^{eco})**

Applications

- ▶ Cell culture aspiration
- ▶ Laboratory waste suction

Ordering Information

197205-01(02)

Lafil 200 - Plus, Suction System with AC100-240V adaptor, US plug (EU plug)

197203-01(02)

Lafil 200 - BioDolphin, Suction System with AC100-240V adaptor, US plug (EU plug)

197215-01(02)

Lafil 200^{eco}- Plus, Suction System with AC100-240V adaptor, US plug (EU plug)

197213-01(02)

Lafil 200^{eco}- BioDolphin, Suction System with AC100-240V adaptor, US plug (EU plug)

197200-53

PP Vacuum Bottle with Quick-Coupling Connector, 4000 mL

197000-60-3500

Laboratory Glass Bottle with GL45 Suction Cap, 3500 mL



* The nameplate is different on Lafil 200 and Lafil 200^{eco}.

Specification

Model	Lafil 200 - Plus / BioDolphin	Lafil 200 ^{eco} - Plus / BioDolphin
PERFORMANCE DATA		
Vacuum range	1013~213 mbar abs.	880~370 mbar abs.
Max. flow rate	10 L/min	
Aspiration rate	17 mL/sec	
Waste bottle	4L PP vacuum bottle / (optional) 3.5L laboratory glass bottle	
Auto shut-off	-	Yes
Noise level	55 dB	
Net weight	2.1 Kg	2.3 Kg
Dimension (LxWxH)	38 x 26 x 22 cm	
ELECTRICAL DATA		
Power input	DC12V, 1.5A	
Max. power	18W	
OTHER DATA		
	Plus	BioDolphin
Suction Kit	<ul style="list-style-type: none"> • Handle • 1-channel tip adaptor (150 mm) • 1-channel tip adaptor with ejector • 8-channel tip adaptor with ejector 	<ul style="list-style-type: none"> • Handle • 1-channel tip adaptor (50 mm) • 1-channel tip adaptor (150 mm) • 1-channel SS adaptor (40 mm) • 1-channel SS adaptor (80 mm) • 8-channel SS adaptor (40 mm) • 1-channel tip adaptor with ejector • 8-channel tip adaptor with ejector • Rack