

Colony Counter

Model: CC-01

Features :

- Contact manual colony counter: a pen gives three sigalsaudible,mark on the plate and addition on the display,anytime that a colony is counted.
- > A model accepts a stainless steel pen tip.
- > It counts by touching the agar directly.
- It can be used for any media just changing the drak field by a transparent grid.



Features :

Readoutz:	4 digit(0-9999)LED display
Mangnification:	2-3X
Petri-dish size:	60-150mm diameter
Zero reset:	yes
Count back:	yes
Dimension:	40x30x30cm