
## SCO -Tech







**REV-100** LCD Digital Rotary Evaporator

## **LCD** display



## **Specifications**

77.1.	100 100 200 010
Voltage [VAC]	100-120/200-240
Frequency [Hz]	50/60
Power [W]	1400
Motor Principle	DC brushless motor
Speed Range [rpm]	20-280
Speed Display	LCD
Reversible Direction of Rotation	Yes
Heating Temperature Range [°C]	RT-180
Heat Control Accuracy [°C]	±1
Heat Output [W]	1300
Temperature Display	LCD
Lift	Motor
Stroke [mm]	150
Timer	Yes
Timer Display	LCD
Time Setting Range [min]	1-999
Overall Dimensions [ D × W × H mm]	465 × 457 × 583
Weight [kg]	15
Permissible Ambient Temperature [°C]	5-40
Permissible Relative Humidity	80%
Protection Class acc. To DIN EN60529	IP20
USB	Yes