


SCODisp

Bottle-top Dispenser



Features:


- Excellent chemical resistance
- Fully autoclavable at 121°C
- Four ranges of bottle-top dispenser cover a volume range from 0.5ml to 50ml
- Easy to clean and maintain
- The optional flexible discharge tube with safety handle permits fast and precise dispensing
- Made of PTFE, FEP, BSG, PP
- Vapor pressure Max. 500mbar, viscosity max. 500mm²/s, temperature max. 40°C, density max. 2.2g/cm³
- SCODisp is supplied with S40, GL32, GL38, GL25, GL28 sized adapters

Specifications

Volume Range [ml]	Graduation [ml]	A ≤ ±		CV ≤ ±	
		%	μl	%	μl
0.5 - 5	0.1	0.5	25	0.1	5
1.0 - 10	0.2	0.5	50	0.1	10
2.5 - 25	0.5	0.5	125	0.1	25
5.0 - 50	1.0	0.5	250	0.1	50

Order Information:

Order No.	Model	Description	
0.382.050	SCOal/sp	Bottle Top Dispenser, 0.5-5ml, incl. adapters	1 ea
0.382.051	SCOal/sp	Bottle Top Dispenser, 1-10ml, incl. adapters	1 ea
0.382.052	SCOal/sp	Bottle Top Dispenser, 2.5-25ml, incl. adapters	1 ea
0.382.053	SCOal/sp	Bottle Top Dispenser, 5-50ml, incl. adapters	1 ea
Accessories:			
0.382.055	S40	Adapter, 45/40 mm	1 ea
0.382.056	GL32	Adapter, 45/32 mm	1 ea
0.382.057	GL38	Adapter, 45/38 mm	1 ea
0.382.058	GL25	Adapter, 32/25 mm	1 ea
0.382.059	GL28	Adapter, 32/28 mm	1 ea
0.382.054	RB	Reagent Bottle (Brown, 1L)	1 ea
0.382.060	FT	Filling tube (m)	1 ea
0.382.061	DT	Discharge tube (m)	1 ea



Adapters



Address. S C O - T e c h
Auf der Heide 15
D-37351 Dingelstädt
Federal Republic of Germany

Phon: +49 (0) 36075 / 43930-6
Fax: +49 (0) 36075 / 43930-8
E-mail: sco@sco-tech.com
Internet: www.sco-tech.com

Contact