

Laboratory Precision Balance

TD- serie

Digital lab Balance, Round or Square pan, Aluminum alloy Shell, LCD Display,
Weighing capacity up to 4000g, Readability: 0.01g, Repeatability: $\pm 0.02g$
For Laboratory or jewelry use

Features:

- High Precision load cell
- Aluminum Alloy die cast base and mains
- Stainless Steel Platter
- High contrast easy to read LED display
- Full Capacity Subtraction
- Simple 3 key panel for easy operation
- Multi Weighing unit conversation g/mg/ct/oz.
- Check Weighting, Piece counting function
- Below Balance Weighing Facility
- Stabilization time 1.5 seconds
- Overload Protection
- RS232 Interface.
- CE approval



Specifications:

Readability:	0.01 g
Repeatability:	$\pm 0.02g$
Linearity:	$\pm 0.02g$
Pan Size	$\varnothing 125mm$ or $170 \times 170mm$
Output Interface	RS232
Dimensions (L x W x H)	$29 \times 18.5 \times 8$ cm
Power Supply	110V/220-240V AC

Model	Weighing Range (g)	Readability Repeatability (g)	Pan Size (mm)	
TD100	100	0.01 $\pm 0.02g$	$\varnothing 125$	
TD200	200			
TD300	300			
TD500	500			
TD600	600			
TD1000	1000			
TD2000	2000		170×170	
TD3000	3000			
TD4000	4000			
TD5000	5000			

SCO Engineering GmbH

Auf der Heide 15
D- 37351 Dingelstädt

Germany

Ph ++49 (36075) 545-70

Fx ++49 (36075) 545-77

E-mail: sco@sco-engineering.com

Internet www.sco-engineering.com