

Type: C-mini

Mini Centrifuge



C-mini Centrifuges performance max. speed of 7000rpm. Ideal for quick spin downs of micro-tubes and PCR tubes, available in four colours green, blue, yellow and pink.

Features

- Supplied with an 8 place 1.5/2.0ml microtube rotor and PCR rotor
- Compact design
- Whisper quiet operation
- Simply close the lid for acceleration and close the lid for deceleration, or use the on/off switch
- Dual door interlock design, for safety
- Multi voltage operation with stable running and high accuracy of speed
- DC motor drive
- Quiet running, noise level <45 dB
- The rotor is seated on a tapered shaft which facilitates easy removal
- Ideal for quick spin downs and micro filtration
- Accelerates and brakes in seconds
- The centrifuge starts and stops by closing and opening lid or using the on/off switch
- Precise control of speed and run time with efficient Separation

Specifications

Model	C-mini
Max Speed	7000 rpm
Max RCF	2680 x g
Rotor Capacity	8×0.2/0.5/1.5/2.0ml tube; 16×0.2ml PCR strips 2×0.2ml PCR strips
	Slide rotor
Run Time	Continuous Operation
Temperature range	ambient
Driving Motor	DC motor
Power	AC100-240 V, 50/60 Hz
Noise level(dB)	≤45
Dimensions	150×150×117 mm [D×W×H]
Weight	0.5 kg



Accessories:

		
Rotor P/N	R8-2	R2-PCR8
Max. speed (rpm)	7000	7000
Max. RCF (× g)	2680	2680
Rotor capacity	2ml/1.5ml x8 0.2ml x 8 0.5ml x 8	16x0.2ml PCR tube 2 × 0.2ml PCR8 tube
Rotor material	High strength plastic	High strength plastic

Order Information:

Item No.	Type	Description
0.154.001	C-mini (b)	Mini Centrifuge with blue lid, incl. R8-2 & R-PCR8
0.154.002	C-mini (g)	Mini Centrifuge with green lid, incl. R8-2 & R-PCR8
0.154.003	C-mini (p)	Mini Centrifuge with pink lid, incl. R8-2 & R-PCR8
0.154.004	C-mini (y)	Mini Centrifuge with yellow lid, incl. R8-2 & R-PCR8
0.154.005	R8-2	Plastic rotor, 7000 rpm, 8 x 2ml/1.5ml
0.154.006	R2-PCR8	Plastic rotor, 7000 rpm, PCR8 x 2
0.154.007	Ad2-02	2ml to 0.2ml rotor adapter, 8 pcs/pk [for R8-2]
0.154.008	Ad2-05	2ml to 0.5ml rotor adapter, 8 pcs/pk [for R8-2]



Address. S C O - T E C H
Auf der Heide 15
D-37351 Dingelstädt
Federal Republic of Germany

Phon: +49 (0) 36075 / 43439-0
Fax: +49 (0) 36075 / 43439-9
E-mail: sco@sco-tech.com
Internet: www.sco-tech.com

Contact