# IECHNOLOGY FOR FUTURE



### Type: CM-mini

# High Speed Mini Centrifuge







Compact, light weight and powerful centrifuge for small sample volumes, complete with rotor 12x 2.0/1.5 ml. It is widely used in physical and chemistry analysis, biochemistry, cellular and molecular biology, clinical Labs and blood donation centers, and is ideal for high-end research applications: Sedimentation of cells and viruses

Separation of sub-cellular organelles

Isolation of macromolecules such as DNA, RNA, proteins or lipids





CM-mini plus

High strength plastic rotor Low noise ≤54dB Delivers fast acceleration and deceleration

#### Features

Powerful centrifugal force 15100 g, wide range of application

CM-mini

- Cast iron chassis with low gravity, provides stable and reliable running at maximum speed, low noise
- User-friendly with large LCD display
- > Brushless motor drives quickly and effortlessly accelerates the rotor to set speed
- > Unique cooling design, minimizes rotor heating
- Fully autoclaveable high strength aluminium alloy rotor
- > Compact design, save valuable space
- > Conform to international safety standards and regulations



## Specifications

CM-mini / CM-mini plus
15.000 rpm (500-15000 rpm), increment: 100 rpm
15100 x g, increment: 100 x g
$12 \times 2.0/1.5$ ml ( or $0.2$ / $0.5$ ml with adapter )
30sec 99 min. – HOLD ( Continuous Operation )
Brushless DC motor
Door interlock. Over-speed detector Over-Temperature detector;
AC 110-220V, 50-60Hz, 1 ph / 3 A
≤ 56 dB [CM-mini] ≤ 54 dB [CM-mini plus]
25s↑ 25s↓ [ CM-mini ] 11s↑ 9s↓ [ CM-mini plus ]
255 x 245 x 140 mm [LxDxH]
6 kg
Speed / Acceleration switch Short-time run function Sound-alert function.



#### **Accessories:**

	Co	
Rotor P/N	R12-2	R12-2P
Max. speed	15000 rpm	15000 rpm
Max. RCF	15100 x g	15100 x g
Rotor capacity	2ml/1.5ml x 12	2ml/1.5ml x 12
	0.2ml x 12	0.2ml x 12
	0.5ml x 12	0.5ml x 12
Rotor material	Aluminium alloy	High strength plastic

#### **Order Information:**

Item No.	Туре	Desctription
0.154.010	CM-mini	High Speed Mini Centrifuge with aluminium rotor R12-2
0.154.011	CM-mini plus	High Speed Mini Centrifuge with plastic rotor R12-2P
0.154.015	R12-2	Aluminium rotor kit, 15000rpm, 12 tubes x 2ml
0.154.016	R12-2P	Plastic rotor kit,15000rpm, 12 tubes x 2ml
0.154.017	Ad2-02	0.2ml rotor adapter, 24pcs/pk
0.154.018	Ad2-05	0.5ml rotor adapter, 24pcs/pk



Federal Republic of Germany

Address. S C O - T e c h
Auf der Heide 15
D-37351 Dingelstädt

Phon: +49 (0) 36075 / 434-390
Fax: +49 (0) 36075 / 434-399
E-mail: sco@sco-tech.com Internet: www.sco-tech.com

Contact